## **Technowood**<sub>®</sub>

HEART LUNG SYSTEM

# Compact Heater-Cooler

AIN

With its compact design, CP4000H can be used for cardiac bypass surgery and extracorporeal membrane oxygenation (ECMO).

#### **FEATURES**


- Temperature in the tank and operation of the water flow pump can be easily operated with the detachable controller. The controller also features an LCD touchpanel for improved visibility and operability. Water circulation state can be confirmed with a single glance from the colors on the screen.
- Although CP4000H is a single-tank heater-cooler device, in heating mode, since only water in the inflow line and the water pump is heated, it is possible for water at 1°C to rise to 36°C within 5 minutes.
- For improved safety and prevention of dry heating, water flow sensors are installed inside the tank and on the inflow line. Heating and cooling will stop if no water flow is detected.



In-house test performance with 3.5L of water. Warming is done through a closed circuit system.

HEATER-COOLER CP4000H	
Temperature setting	0-42°C (resolution: 0.1°C)
Туре	Cooler: 200WCompressor Heater: Built-in 1kW
Water flow pump	Max: 11L/min
Power supply	AC100V - 230V; 50/60Hz
Power consumption	1100VA
Tank volume	2.5L
Dimension (WxDxH) *without mast	280mm×550mm×517mm 38kg

#### **DEVICE PERFORMANCE**






The LCD touchpanel can be detached from the unit and attached to a mast closer to the surgical field.



### CP4000H (T)

Models equipped with panels and other accessories to accomodate portable extracorporeal membrane oxygenation (ECMO) systems are also available. **CP4000H HEATER-COOLER ACCESSORIES** 



#### JAPAN

#### **Technowood Corporation**

Tel: +81 (3) 3856-4111 Fax: +81 (3) 3856-4113

#### USA

#### **Technowood America Corporation**

Tel: +1 (714) 434-8713 Fax: +1 (714) 434-8715

#### INTERNATIONAL

#### **Technowood International Pte. Ltd.**

Tel: +81 (3) 3898-5252 Fax: +81 (3) 3898-5252

Please refer to the Instructions For Use for details. Terms and conditions may change without prior notice. Distributed by:



www.technowood.co.jp

BS1912R3